

Infectious Disease Collection Guide – Validated Sample Types



UTM Swab (MicroTest M6 Kit)

Micro-tipped and Traditional-Tipped swabs

- Respiratory Assays
 - COVID 19
 - RPP panels
 - Bordetella
- HSV 1+2/VZV



UTM Swab

- Women’s Health (also above)
 - Bacterial Vaginosis - I
 - Bacterial Vaginosis - II
 - Candidiasis Panel
 - STD Panel
 - *T. vaginalis*
 - *M. genitalium*



Aptima Multitest Swab Collection Kit

- STI Panel (females only)
 - *C. trachomatis/N. gonorrhoeae*
 - *T. vaginalis*
- *M. genitalium* (males and females)
- Women’s Health (also below)
 - Aptima BV Assay
 - Aptima CV/TV Assay



Urine with Boric Acid

- Urine Culture



Aptima Urine Specimen Collection Kit (males and females)

- STI Panel *C. trachomatis/N. gonorrhoeae*
 - *T. vaginalis*
- *M. genitalium*

Note: urine must be transferred from a sterile container into the collection kit within 24 hours



Orange Cap Fecal Swab

- Gastrointestinal Infection



Throat Swab (double swab)

- Strep Complete